

MEGADYNE

Generator Compatibility Chart Updated 6/21/04

Compatible with the MEGA 2000 Patient Return Electrode System

Manufacturer	Model	Isolated - (Most are Compatible w/Mega 2000 see comments)	Adapter or Detachacable Required for MEGA 2000 other than M2K- 01	Comments
Aspen	MF380	Yes	M2K-03, M2K-04	
Aspen	MF360A	Yes	M2K-03, M2K-04	
Aspen	MF450	Yes	M2K-03, M2K-04	
Aspen/Conmed	Sabre 2400	Yes	M2K-03, M2K-04	
Aspen/Conmed	Sabre180	Yes	M2K-03, M2K-04	
Aspen/Conmed	Excalibur Plus	Yes	No	
Bard	3000	Yes	M2K-03, M2K-04	
Berchtold	Elektrotom 400	Yes	M2K-03, M2K-04	
Berchtold	Elektrotom 640	Yes	M2K-03, M2K-04	
Berchtold	Elektrotom 200	Yes	M2K-03, M2K-04	
Berchtold	Elektrotom 200-165	Yes	M2K-03, M2K-05	
Berchtold	Elektrotom 300B	Yes	M2K-03, M2K-04	
Berchtold	Elektrotom 390	Yes	M2K-03, M2K-04	
Berchtold	Elektrotom 621	Yes	M2K-03, M2K-04	
Birtcher	6400	Yes	A-238, M2K-06	Argon Beam Coagulator
Birtcher	4400	Yes	No	
Birtcher	6000	Yes	See comments	Argon Beam Coagulator - this model was originally designed for single plate return electrodes and will work with Mega 2000. Some units have been retrofitted with PSS (pad sensing system) and do not work with Mega 2000.
Birtcher/Bard/EMS/Davol/NDM	5000	Yes	M2K-03, M2K-04	
Circon	BC-200	Yes	Yes	
Circon	91-I	Yes	Yes	
Concept	9900	Yes	M2K-03, M2K-04	
Concept	9700	Yes	M2K-03, M2K-04	
Concept	9600	Yes	M2K-03, M2K-04	
Conmed	6500	Yes	A-238, M2K-06	Argon Beam Coagulator
Conmed	7500	Yes	A-238, M2K-06	Argon Beam Coagulator
ConMed	System 5000	Yes	M2K-01 or M2K-06	The shape of the return electrode receptacle may require the M2K-06 be used.
Cooper	Leep 1000	Yes	M2K-03, M2K-04	
Erbe	ICC350	Yes	No	Do not use in with the High Cut or Endo Cut Mode. Doing so may result in a greater electrosurgical effect than intended.
Erbe	ICC300	Yes	No	Do not use in with the High Cut or Endo Cut Mode. Doing so may result in a greater electrosurgical effect than intended.
Erbe	ICC200	Yes	No	Do not use in with the High Cut or Endo Cut Mode. Doing so may result in a greater electrosurgical effect than intended.
Eschmann	TD 311	Yes	M2K-03, M2K-04	
Eschmann	TD 302	Yes	M2K-03, M2K-04	
Eschmann	TD 411	Yes	M2K-03, M2K-04	
Eschmann	TD 402	Yes	M2K-03, M2K-04	

Manufacturer	Model	Isolated - (Most are Compatible w/Mega 2000 see comments)	Adapter or Detachacable Required for MEGA 2000 other than M2K- 01	Comments
Eschmann	TD 411RS	Yes	M2K-03, M2K-04	
Eschmann	TD411RS2	Yes	M2K-03, M2K-04	
Eschmann	TD 830	Yes	M2K-03, M2K-04	
Eschmann	TD 850	Yes	M2K-03, M2K-04	
Ethicon (Formerly from Medtrex)	Pegasys (OR Pro)	Yes	No	This generator has been discontinued and should be removed from the field.
Everest	8750	Yes	No	? With MEGA 2000 special procedure unit, has not been tested with MEGA 2000
Martin	ME MB1	Yes	M2K-03, M2K-04	
Martin	ME 200	Yes	M2K-03, M2K-04	
Martin	ME 411	Yes	M2K-03, M2K-04	
Martin	ME-300	Yes	M2K-05	
Martin	ME-400	Yes	M2K-05	
Maxxim/Bovie	X-15	Yes	No	
Maxxim/Bovie	X-30	Yes	No	
Maxxim/Bovie	X-40	Yes	No	
NDM	1000 Powerpoint	Yes	M2K-03, M2K-04	
Neomed	300A	Yes	M2K-03, M2K-04	Unit was made by Birtcher
Olympus	UES	Yes	Yes	Must obtain from Olympus
Olympus	PSD-20	Yes	Yes	Must obtain from Olympus
Orotec Interventions	Ora 50	Yes	Yes, adapter 805019 from Orotec	Dedicated to Arthroscopy - Call MegaDyne re: questions
Richard Wolf	2083	Yes	???	
Richard Wolf	2093	Yes	???	
Richard Wolf	2094	Yes	???	
Richard Wolf	2353	Yes	???	
Richard Wolf	4083	Yes	???	
Valleylab	Surgistat	Yes	E-0504-1L	Adapter required on older units only
Valleylab	SSE2L	Yes	E-0504-1L	
Valleylab	Force FX	Yes	No	
Valleylab	Force 2	Yes	No	
Valleylab	Force 2 CEM	Yes	No	
Valleylab	Force 1C	Yes	No	
Valleylab	Force 1B	Yes	No	
Valleylab	Force 30	Yes	No	
Valleylab	Force 40	Yes	No	
Valleylab	Force EZ	Yes	No	

Manufacturer	Model	Isolated - (Most are Compatible w/Mega 2000 see comments)	Adapter or Detachacable Required for MEGA 2000 other than M2K- 01	Comments
NOT Compatible with the MEGA 2000 Patient Return Electrode System				
Orotec Interventions	Vulcan	Isolated - not compatible with Mega 2000	No	Dedicated to Arthroscopy - functions with monitoring pads only
Alcon	203-0000-501	NA		Ultrasonic Device - does not use a return electrode
Codman	80-1170	NA		Bipolar Only Machine - No Return Electrode Required
Ethicon	G110	NA		Ultrasonic Scalpel - does not use a return electrode
Valleyforge	Symmetry	NA		Bipolar Only Machine - No Return Electrode Required
Maxxim/Bovie	400SR	No	????	Not compat. w/MEGA 2000
Maxxim/Bovie	CSV Bovie	No	B-205	Not compat. w/MEGA 2000
Maxxim/Bovie	URO Bovie	No	B-205	Not compat. w/MEGA 2000
Maxxim/Bovie	Bantam/Bovie	No	E-0504-1L	Not compat. w/MEGA 2000
Maxxim/Bovie	Ritter A	No	E-0504-1L	Not compat. w/MEGA 2000
Maxxim/Bovie	X-10	No	No	Not compat. w/MEGA 2000
Maxxim/Bovie	Specialist	No		Not compat. w/MEGA 2000
Valleylab	SSE3	No	E-0504-1L	Not compat. w/MEGA 2000
Valleylab	SSE3B	No	E-0504-1L	Not compat. w/MEGA 2000
Valleylab	Force 4	No	No	Not compat. w/MEGA 2000
Valleylab	Force 4B	No	No	Not compat. w/MEGA 2000
Cryomedics	91-J	Yes	???	? With MEGA 2000 special procedure unit, has not been tested with MEGA 2000
EP Technologies	EPT-1000	Yes	No	? With MEGA 2000 special procedure unit, has not been tested with MEGA 2000
Medtronic	601	Yes	???	? With MEGA 2000 special procedure unit,
Medtronic	Cardorhythm	Yes	???	? With MEGA 2000 special procedure unit, has not been tested with MEGA 2000
Microvasive	Endostat	Yes	E-0504-2	? With MEGA 2000 special procedure unit, has not been tested with MEGA 2000
Radionics	RFG-3C	Yes	???	? With MEGA 2000 special procedure unit, has not been tested with MEGA 2000
Ellman	Surgitron			No adapters available for this machine