
Certificate of Compatibility

Megadyne Medical Products Inc. herewith declares, that our product, manufactured under the requirements of 93/42/EEC in combination with EN 60601 (Certificate Q252.552), and the registered Quality System Registration NSAI (0050) MD19.2573:

0800	MEGA 2000® Reusable Patient Return Electrode
0830	MEGA Soft® Patient Return Electrode
0835	MEGA Soft Dual Cord Patient Return Electrode
0840	MEGA Soft Pediatric Patient Return Electrode

In combination with

M2K-08	MEGA Soft Reusable Compatibility Cable (8.0 ft / 2.4 m)
--------	--

OR

M2K-09	MEGA Soft Reusable Compatibility Cable (14.5 ft / 4.4 m)
--------	---

is compatible with the electrosurgical generator

Manufacturer: **Valleylab**

Model Code : **Force Triad™**

Draper, Utah, November 2009



Ronda Magneson
Director, Regulatory Affairs
Megadyne Medical Products